



# THE TRUE SHAPE OF IOP

The CATS® Tonometer Prism is designed to establish a new standard in reliable Goldmann IOP measurement accuracy by addressing key corneal biometrics including Center Cornea Thickness (CCT), Corneal Hysteresis (CH), Corneal Curvature and Tear Film which affect 50% of the population. Accurate IOP measurements provide eye care clinicians with a key metric in the diagnosis and treatment of glaucoma and other serious ocular conditions.

Our technology provides better information and instills confidence in the IOP measurements used to assess and protect your patient's sight.


The CATS® Tonometer Prism's advanced technology coupled with its improved ergonomic function and aesthetic appeal provides complete clinical ease of use. We believe the CATS® Tonometer Prism will quickly supplant traditional prisms to become the new standard of care for IOP measurement for all eye care clinicians (ophthalmology and optometry) and all patient populations.

How often do you ponder these questions (do these questions trouble you)?

**“Does my patient have glaucoma?”**  
**“Is the pressure low enough?”**

Available in both reusable and disposable options, the FDA cleared and patented CATS® Tonometer Prism is designed to deliver improved accuracy for a majority of the population. The CATS® Prism will seamlessly integrate into a clinician's practice with no additional cost and no new training...while significantly streamlining patient care.

## What Others are saying about the CATS® Prism



"The CATS® Tonometer Prism is designed to provide more accurate IOP assessment, sparing patients and clinicians from the damaging consequences of over and under estimation of IOP."

- Nathan Radcliffe, MD



"The CATS® Tonometer Prism is designed to deliver more accurate IOPs, because it compensates for corneal biomechanics."

- John Berdahl, MD

Available in Reusable and Disposable, Single-Use options:



**Reusable**

**Item number: CATS-R**

Ergonomic design for exceptional ease of use.  
Designed for improved accuracy.  
Fits all Goldmann and Perkins tonometers.



**Disposable, Single-Use (100 prisms/box)**

**Item Number: CATS-D**

Featuring Snap-Fit™ System for one-handed attachment and disposition of tonometer tips.  
Ergonomic Design for exceptional ease of use.  
Designed for improved accuracy.  
Fits all Goldmann and Perkins tonometers.

## AWARDS



CAT04106 Rev A



## CONTACT US

1517 N Wilmot Rd, Suite 143, Tucson, AZ 85712  
customerservice@catsiop.com | 520-370-4113 | catsiop.com



Virtual Press Kit:

[catsiop.com/info](https://catsiop.com/info)